



# K&S Associates, Inc.

1926 Elm Tree Drive  
Nashville, Tennessee 37210-3718  
800-522-2325 Fax 615-871-0856



## A2LA ACCREDITED INSTRUMENT CALIBRATION LABORATORY

FACILITY \_\_\_\_\_ PHONE \_\_\_\_\_

PERSON \_\_\_\_\_ FAX \_\_\_\_\_

INSTRUMENT \_\_\_\_\_

### A2LA ACCREDITED PRECISION ION CHAMBER CALIBRATION:

\*\*\*\*\*

Energy	Range	Uncertainty
Cobalt 60 air kerma	30-1000 mGy/min	1.5%
Cesium 137 air kerma	5-120 mGy/min	1.6%
M series X-rays	0.3 – 200 mGy/min	1.2%
M30 _____	M100 _____	
M40 _____	M120 _____	
M50 _____	M150 _____	
M60 _____	M175 _____	
M70 _____	M200 _____	
M80 _____	M220 _____	
M90 _____	M250 _____	
NS Series X-Rays	6-1000 uGy/min	1.2%
NS40 _____	NS150 _____	
NS60 _____	NS250 _____	
NS80 _____		
NS100 _____		

### A2LA ACCREDITED ENVIRONMENTAL LEVEL ION CHAMBERS:

\*\*\*\*\*

Energy	Range	Uncertainty
Cesium 137	0.4 – 5 uGy/h	11%

### A2LA ACCREDITED ELECTROMETER CALIBRATIONS:

\*\*\*\*\*

Energy	Range	Uncertainty
Charge &/or Current	any	0.1%

**Calibration intervals: Re-calibration reminder service will be set on a 12 month interval unless you specify otherwise. If you wish a re-calibration date recorded on your report please indicate the interval you require in the space provided.**

American Express, MasterCard, Visa are accepted email: ksinfo@kslab.com web: http://www.kslab.com