



K&S Associates, Inc.
 1926 Elm Tree Drive
 Nashville, Tennessee 37210-3718
 800-522-2325 Fax 615-871-0856



HPS ACCREDITED INSTRUMENT CALIBRATION LABORATORY

FACILITY _____ PHONE _____

PERSON _____ FAX _____

INSTRUMENT _____

HPS ACCREDITED PRECISION ION CHAMBER CALIBRATION:

| Energy | Range | Uncertainty |
|----------------------|-------------------|-------------|
| Cobalt 60 air kerma | 30-1000 mGy/min | 1.5% |
| Cesium 137 air kerma | 5-120 mGy/min | 1.6% |
| M series X-rays | 0.3 – 200 mGy/min | 1.2% |
| | M30 _____ | |
| | M50 _____ | |
| | M80 _____ | |
| | M100 _____ | |
| | M150 _____ | |
| | M200 _____ | |
| | M250 _____ | |
| NS Series X-Rays | 6-1000 uGy/min | 1.2% |
| | NS80 _____ | |
| | NS100 _____ | |
| | NS150 _____ | |
| | NS200 _____ | |

HPS ACCREDITED ENVIRONMENTAL LEVEL ION CHAMBERS:

| Energy | Range | Uncertainty |
|------------|---------------|-------------|
| Cesium 137 | 0.4 – 5 uGy/h | 11% |

A2LA ACCREDITED ELECTROMETER CALIBRATIONS:

| Energy | Range | Uncertainty |
|---------------------|-------|-------------|
| Charge &/or Current | any | 0.1% |

Calibration intervals: Re-calibration reminder service will be set on a 12 month interval unless you specify otherwise. If you wish a re-calibration date recorded on your report please indicate the interval you require in the space provided.

American Express, MasterCard, Visa are accepted email: ksinfo@kslab.com web: <http://www.kslab.com>