



**K&S Associates, Inc.**  
 1926 Elm Tree Drive  
 Nashville, Tennessee 37210-3718  
 800-522-2325 Fax 615-871-0856



**KEITHLEY 35050A/35080 SYSTEM**

FACILITY: \_\_\_\_\_

CONTACT PERSON: \_\_\_\_\_

FAX: \_\_\_\_\_ PHONE: \_\_\_\_\_

**SPECIAL Kit Calibration**  
 includes the following:

**precalibration evaluation of all kit components**

**35080 kVp divider** noninvasive kVp meter calibration, includes all filter packs

**diagnostic range** - W target, Al filter, 50 - 150 kVp, (3.48 mm Al HVL at 80 kVp)

**mammography range** - Mo target, Mo filter, 22 - 35 kVp ( .35mm Al @ 28 kVp)

**35050 dosimeter** calibration includes entry of calibration factors for each chamber

**96035 15cc chamber ADCL calibration**  
 28 kVp Mo/Mo (MQSA compliant)  
 100 kVp general rad

**96020 150cc chamber ADCL calibration**  
 100 kVp general rad

timer accuracy

mA/mAs meter accuracy

**Calibration intervals:**

Your re-calibration reminder will be set up on a "manufacturer's recommended" interval unless you specify otherwise. IF you wish a re-calibration date recorded on your report please indicate the interval you require in the space below.

American Express, Visa or MasterCard are accepted.  
 Shipping and handling is not included in the above fees.

Phone 800-522-2325 Fax 615-871-0856  
 email: [ksinfo@kslab.com](mailto:ksinfo@kslab.com) web: <http://www.kslab.com>