



K&S Associates, Inc.

1926 Elm Tree Drive
Nashville, Tennessee 37210-3718
800-522-2325 Fax 615-871-0856



NON-INVASIVE kVp METER **A2LA ACCREDITED CALIBRATION**

Facility: _____ Phone: _____

Person: _____ Email: _____

Instrument Description: _____

Calibration Units and Available Ranges:

Diagnostic Radiology: W target, Al filter x-ray tube, 50 – 150 kVp, 3.4 mm Al @ 80 kVp

Mammography: Mo target, Mo filter x-ray tube, 22 – 35 kVp, 32 μm Mo filter

Ordering: Please fill out one form for each kVp meter to be calibrated.

Standard Calibration includes 6 individual kVp calibration points across on range.

Timer Accuracy _____

mA/mAs meter Accuracy _____

Dose (relative dose) Linearity _____

External Chamber Calibration (First Point) _____

Additional Calibration Point needed _____

Calibration intervals: Your re-calibration reminder will be set on a 24 month interval unless you specify otherwise. IF you wish a re-calibration date recorded on your report please indicate the interval you require in the space below.

American Express, MasterCard, Visa accepted

email: ksinfo@kslab.com

web: <http://www.kslab.com>