



**K&S Associates, Inc.**  
 1926 Elm Tree Drive  
 Nashville, Tennessee 37210-3718  
 800-522-2325 Fax 615-871-0856



**BRACHYTHERAPY  
 ACCREDITED DOSIMETRY CALIBRATION LABORATORY**

Facility \_\_\_\_\_ Phone \_\_\_\_\_

Person \_\_\_\_\_ Fax \_\_\_\_\_

WELL CHAMBER CALIBRATION

HIGH DOSE RATE POINT      \_\_\_\_\_ Iridium 192

LOW DOSE RATE POINTS

Iodine 125: -----

- \_\_\_\_\_ Nycomed Amersham Model 6711
- \_\_\_\_\_ Nycomed Amersham Model 6733 Echoseed
- \_\_\_\_\_ Best Ind. Model 2301
- \_\_\_\_\_ Draximage LS-1 (sold by Cytogen)
- \_\_\_\_\_ Ibt 1251
- \_\_\_\_\_ Imagyn Isostar Model IS 12501 (bought by Bard Brachytherapy)
- \_\_\_\_\_ Implant Science Corp 3500
- \_\_\_\_\_ IsoAid Advantage
- \_\_\_\_\_ Mills Biopharmaceuticals Inc. Model I 125SH Prostaseed (now Core MBI)
- \_\_\_\_\_ BSI Brachytherapy IOGOLD Model MED 3631 A/M
- \_\_\_\_\_ SourceTech Medical STM 1251 (now Bard Brachytherapy)
- \_\_\_\_\_ Bebig (Isotope Products) Model I125.S06

Palladium103: -----

- \_\_\_\_\_ Best Ind, Model 2335
- \_\_\_\_\_ Draximage Pd-1
- \_\_\_\_\_ Int'l Brachytherapy Model 1031L
- \_\_\_\_\_ BSI Brachytherapy Model MED 3633 A/H
- \_\_\_\_\_ Theragenics Theraseed Model MO 200
- \_\_\_\_\_ Bebig (Theragenics) Model Pd 3.S11 (sold by Theragenics)
- \_\_\_\_\_ Mills

Cesium 131:-----

- \_\_\_\_\_ IsoRay Cs-1

Iridium 192: -----

- \_\_\_\_\_ Best Ind. Iridium 192 85-01

Cesium 137: -----

- \_\_\_\_\_ 3M Cesium 137

**Calibration intervals: Your re-calibration reminder will be set on a 24 month interval unless you specify otherwise. IF you wish a re-calibration date recorded on your report please indicate the interval you require in the space below.**

BRACHYTHERAPY SOURCE(S) CALIBRATION \_\_\_\_\_  
 Shipping container rental \_\_\_\_\_

email: [ksinfo@kslab.com](mailto:ksinfo@kslab.com)    web: <http://www.kslab.com>